

Body bag with window




Close



Open

【Specification】

Catalogue #	Y-3001
Size	W : 800 × L : 250mm
Material	Polypropylene, 140g / m ²
Colour	White
Load capacity	160kg
Package	1 sheet
Delivery unit	10 sheets / Box



Normeca Asia Co., Ltd

19-3, Asakusabashi 3-chome, Taito, Tokyo 111-0053
TEL: +81-3-5687-3899, FAX: +81-3-5687-3911
Email: sales@normeca-asia.jp