

HOSPITAL BED

The NorCare Bed

 **Normeca** AS
SUPPLIER OF MOBILE HOSPITALS
Mallandvn 12, P.O. Box 404, 1470 Lorenskog - Norway
Tlf.: +47 67 92 76 00 - Fax: +47 67 92 76 92
Mail: post@normeca.no - Web: www.normeca.no



THE NORCARE BED

Equally suited for emergency conditions as for long-term care

THE NORCARE BED IS A LIGHT WEIGHT BED, CAPABLE OF ENDURING MAXIMUM STRESS AND WEIGHT. IN A FIELD HOSPITAL SETTING, THE HIGH QUALITY OF THE NORCARE BED MAKES IT SUITABLE FOR INTERNAL MEDICAL- AS WELL AS FOR SURGICAL INPATIENTS.

The add-on equipment allows optimal areas of use. The bed is equally suited for emergency conditions as for more long-term care. The height of 67 cm between floor and patient area facilitates getting in and out of bed for patient, and gives medical staff a good working position.

A reserve supply of the NorCare Bed is also the perfect solution in emergency response units of hospitals, clinics, camps and other accommodation facilities.

The ingenious, yet simple, design allows great flexibility and rational use of storage space.

The NorCare Bed measures only 5 cm in height before assembly. When looking to keep shipping costs down, this is an obvious advantage.

OPTIONAL EQUIPMENT:

- Wheels with brakes
- Traction device
- Night table
- Adjustable back support
- Adjustable infusion stand
- Reading lamp
- Mattress/mattress cover
- Safety railing



QUALITY DESIGN:

- Anodized aluminium is used for all bearing construction.
- The design allows assembly with no previous instructions necessary.
- All moving parts are incorporated in the structure itself, eliminating the risk of losing vital components, such as screws and washers.
- With canvas made from flame-resistant synthetic fabric, the NorCare Bed is hygienic and easy to keep clean.
- Because of the patented aluminium profiles, the accessories are easily attached.



WEIGHT:

16 kg


MEASUREMENTS BEFORE ASSEMBLY:

209 cm x 75 cm x 5 cm



...a light weight bed, capable of enduring maximum stress and weight.

TECHNICAL SPECIFICATIONS:



OPTIONAL EQUIPMENT	WEIGHT	DESCRIPTION	DIMENSIONS
A. Wheels with brakes	2,80 kg	Total length of bed	209 cm
B. Traction device	2,50 kg	Length of patient area	200 cm
C. Night table	1,50 kg	Total breadth of bed	75 cm
D. Adjustable back support	2,50 kg	Breadth inside of legs	56 cm
E. Adjustable infusion stand	1,00 kg	Total height, assembled	93 cm
F. Reading lamp	1,50 kg	Distance, ground to patient	67 cm
G. Mattress/mattress cover	4,50 kg	Area	

MEASUREMENTS OF BED, UNASSEMBLED:

209 cm x 75 cm x 5 cm

WEIGHT OF BED:

16 kg

MEDIQ NORGE:

Preparedness Department is well established in field of medical relief materials to emergency-, disaster-, and developing areas around the world. Mediq Norge has gained wide recognition among national and international NGOs as well as in the UN system and with local authorities as the supplier of Norwegian Mobile Hospital and Disaster Unit (NorHosp) and Norwegian Medical Aid System (NorAid). They are continuously seeking new areas for product development.